

ABC Model 4100

Syringe Infusion Pump



Ease of Operation Particularly for Enteral Feedings

User Customization promotes friendly programming

Multiple Delivery Modes

- Various Delivery modes (Continuous, Mass, Body-Weight, Volume over Time, Bolus option)
- Adjustable Flow Rate ranges
- Accuracy +/- 3% excluding syringe variations

Storage for User Defined Settings

- 90 Settings can be saved with Customized Name

Accommodate Various Syringes Manufacturers with Automatic Size Sensing

- Enteral - 1ml to 140ml
- IV – 1ml to 140 ml

Easy Syringe Loading

- Front Syringe Loading with Syringe Light
- Easy Plunger Loading with Plunger Sensor

Multiple Sensors to Monitor Delivery Accuracy

- Syringe Size
- Syringe Plunger Force
- Syringe Plunger Movement
- Motor Shaft Rotation

Metal Housing Design

- Non-plastic housing to reduce breakage
- Improved EMI immunity (IEC 60601-1-2)
- Splash Proof for more effective and easier cleaning



Large Graphic LCD Display with White Backlight (promotes easy viewing)

Syringe Light can be viewed from afar (delineates infusion status)

Unique Labeling area can be used to indicate infusion information

AC / Battery operated device with Externally Removable Battery Pack

- Standard Battery Capacity (12 hours at 5 ml/hr)
- No External Power Adaptor required

LED & LCD Alerts and Audio Alarms

RS232 Interface Capabilities

Flexible Pole Mount System

Distributed By:



Ph: 847-891-6210 Fax: 847-891-6215