

# WaterPAP®

Bubble PAP Setup  
for Infants and Neonates

## Convenience

Equipment ready  
when needed

## Uniformity

Each patient setup  
exactly the same

## Accuracy

Visual representation of  
prescribed pressure

## Savings

Effective performance  
with lower cost



[www.airwaysdevelopment.com](http://www.airwaysdevelopment.com)

Ordering Information on Reverse Side

Latex Free • Made in the USA • MK-9

©Copyright 2012 Airways Development LLC, All Rights Reserved  
WaterPAP is a Registered Trademark of Airways Development LLC



### Ordering Information

Product #	Description	Quantity
WP-7700	WaterPAP Setup Includes: Water-Tight Canister with 0-10cm Graduated Label, 10mm Corrugated Tube Affixed with Air Diffuser, Tube Lock Top & Protective Cover	Box of 10 Case of 40
WP-7705	WaterPAP Starter Kit, includes: Four WaterPAP Setups, One Pole Clamp, One Bracket with Tube Guide	Each
CP-2203	Pole Clamp	Each
CP-2204	Bracket	Each
WP-7704	Bracket with Tube Guide	Each
CP-2209	22mm ID x 6mm OD Fitting Combination, Connects Oxygen Supply Tube to Humidifier	Box of 9
CP-2215	22x15x6mm TEE Assembly, Accepts 15mm Oxygen Sensor, Supply Tube and Humidifier	Box of 5 Case of 20
CP-2218	22x18x6mm TEE Assembly, Accepts 18mm Oxygen Sensor, Supply Tube and Humidifier	Box of 5 Case of 20

[www.airwaysdevelopment.com](http://www.airwaysdevelopment.com)