





The Plak-Vac® Petite LA™ low absorption swab was developed for applications where the loss of oral therapy solution (e.g., colostrum) must be minimized. The Petite LA swab is size appropriate for even the smallest babies and features a soft, non-shedding foam head that will not lint or fray when wet, and a smooth polypropylene handle rather than wood.



Due to its low absorption and small head, the Petite LA swab is ideal for applications such as the oropharyngeal administration of colostrum.

The Petite LA swab can be incorporated into a Plak-Vac oral care kit to support a protocol for the oropharyngeal administration of colostrum, by combining the swab with other components such as a 1 cc syringe, gauze pads, sterile water, gloves, and other accessories.

Disclaimer: Trademark Medical makes no claim of safety or efficacy for any products used in Plak-Vac Oral Care Kits beyond those made by the product manufacturer's legally approved claims and indications for use.