




Peripherally Inserted Central Catheter with Peel-Away Introducer Sheath



The Footprint™ introducer sheath has unique design features to enhance first stick success.



The hub is clearly marked with "PICC" and French size to readily identify catheter type.

Peripherally Inserted Central Catheter

- Designed especially for neonatal/pediatric patients
- Silicone catheter with integrated extension set
- Soft stabilization wing for enhanced securement
- Indelible markings every centimeter
- Fully radiopaque
- Latex and DEHP free

Introducer Sheath

- Notched needle with minimal lie distance for immediate visualization of flashback designed to increase first stick success
- Low profile, compact introducer designed to enhance small vein access



Footprint™
medical inc.

210.226.2600

www.footprintmed.com



Ordering Information for:

Peripherally Inserted Central Catheters

Silicone Peripherally Inserted Central Catheters

S1PIC1.9-S	Silicone PICC Complete Kit	1.9 Fr.	5/cs
SHEATH1.9	PICC Sheath Introducer Only	1.9 Fr.	5/cs
S1PIC1.9-C	Silicone PICC Catheter Only	1.9 Fr.	5/cs

Catheter					Flow Rates	
Product Number	Size	Usable Length	Gauge	Prime Volume	Pump	Gravity
S1PIC1.9-S	1.9 FR.	28 CM	23 GA.	0.09 ML	240 ML/HR	40 ML/HR

"Helping tiny babies make their own footprints in the world..."

Distributed By:



1203 Camden
San Antonio, Texas 78215
210.226.2600
www.footprintmed.com